

Product datasheet for **RR211262L3V**

Hopx (NM_133621) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hopx (NM_133621) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hopx
Synonyms:	Gllg15b; Hod; Obl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133621
ORF Size:	219 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211262).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133621.2 , NP_598305.2
RefSeq Size:	1065 bp
RefSeq ORF:	222 bp
Locus ID:	171160
UniProt ID:	Q78ZR5
Cytogenetics:	14p11



[View online »](#)

Gene Summary:

Atypical homeodomain protein which does not bind DNA and is required to modulate cardiac growth and development. Acts via its interaction with SRF, thereby modulating the expression of SRF-dependent cardiac-specific genes and cardiac development. Prevents SRF-dependent transcription either by inhibiting SRF binding to DNA or by recruiting histone deacetylase (HDAC) proteins that prevent transcription by SRF. Overexpression causes cardiac hypertrophy. Acts as a co-chaperone for HSPA1A and HSPA1B chaperone proteins and assists in chaperone-mediated protein refolding.[UniProtKB/Swiss-Prot Function]