

Product datasheet for **RR211276L3V**

Cnot10 (NM_001007003) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cnot10 (NM_001007003) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cnot10
Synonyms:	MGC94752
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007003
ORF Size:	2232 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007003.1 , NP_001007004.1
RefSeq Size:	2662 bp
RefSeq ORF:	2235 bp
Locus ID:	316034
UniProt ID:	Q5XIA4
Cytogenetics:	8q32



[View online »](#)

Gene Summary:

Component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Is not required for association of CNOT7 to the CCR4-NOT complex (By similarity).[UniProtKB/Swiss-Prot Function]