

Product datasheet for **RR211278L3V**

Gstk1 (NM_181371) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gstk1 (NM_181371) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gstk1
Synonyms:	GST13-13; GSTkappa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181371
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211278).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181371.2 , NP_852036.1
RefSeq Size:	897 bp
RefSeq ORF:	681 bp
Locus ID:	297029
UniProt ID:	P24473
Cytogenetics:	4q24
Gene Summary:	member of a family of enzymes that catalyse the conjugation of glutathione (GSH) to a variety electrophiles [RGD, Feb 2006]



[View online »](#)