

Product datasheet for **RR211294L4V**

Fam46b (NM_001107907) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam46b (NM_001107907) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fam46b
Synonyms:	RGD1308967
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001107907
ORF Size:	1281 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211294).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001107907.2 , NP_001101377.1
RefSeq Size:	1724 bp
RefSeq ORF:	1284 bp
Locus ID:	313019
UniProt ID:	B0BNK8
Cytogenetics:	5q36
Gene Summary:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. [UniProtKB/Swiss-Prot Function]



[View online »](#)