

## Product datasheet for **RR211338L4V**

### Degs1 (NM\_053323) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Degs1 (NM_053323) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Degs1
Synonyms:	Degs
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053323
ORF Size:	969 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211338).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053323.2</a> , <a href="#">NP_445775.2</a>
RefSeq Size:	1466 bp
RefSeq ORF:	972 bp
Locus ID:	58970
UniProt ID:	<a href="#">Q5XIF5</a>
Cytogenetics:	13q26
Gene Summary:	human homolog is a member of a sphingolipid delta 4-desaturase enzyme family [RGD, Feb 2006]



[View online »](#)