

Product datasheet for **RR211369**

Eif4a1 (NM_199372) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eif4a1 (NM_199372) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Eif4a1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR211369 representing NM_199372 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCGGAGTCAGGATTCTCGATCCAGAGACAATGGCCCCGATGGGATGGAGCCGGAAGCGTCATCG
AGAGTAACGGGAATGAGATTGTGGATAGCTTCGATGACATGAATCTCAGAATCCCTCCCGTGGTAT
TTATGCCTATGGTTTTGAGAAGCCCTCTGCCATCCAGCAGCGAGCTATTCTTCTTGATCAAGGGTTAT
GATGTGATTGCTCAAGCCAGTCTGGGACTGGGAAAACAGCTACATTTGCCATATCCATTCTGCAGCAGA
TTGAATTAGATCTAAAGGCCACTCAGGCTTTGGTTCTGGCACCCACTCGTGAATTGGCTCAGCAGATTCA
AAAGGTGGTTATGGCACTGGGAGACTACATGGGTGCCTCTGTGCATGCCTGTATTGGGGCACCACCGTG
CGTGCTGAGGTGCAGAAGCTGCAAAATGGAAGCTCCCCATATCATCGTGGGCACCCCTGGCCGGGTGTTCC
ATATGCTTAACCGGAGATACCTGTCACCCAAATACATCAAGATGTTGTCCTGGATGAAGCAGATGAAAT
GCTAAGCCGTGGGTTCAAGGATCAGATCTATGATATATCCAAAAGCTCAACAGCAACACACAGGTAGTT
TTGCTGTCTGCCACAATGCCTTCTGATGTCCTTGAGGTGACTAAGAAGTTCATGAGAGACCCTATTCGGA
TTCTTGCAAGAAGGAAGAGTTGACCCTGGAGGGTATCCGCCAATCTACATCAATGTGGAACGAGAGGA
GTGGAAGCTCGACACATTGTGTGACTTGTATGAAACACTGACCATCACCCAGGCAGTAATCTTTATTAAT
ACCAGAAGGAAGGTGGACTGGCTCACTGAGAAGATGCATGCCCGGATTTACAGTTTCTGCCATGCATG
GAGATATGGACAAAAGGAACGAGACGTGATCATGAGGGAGTTCGGTCTGGCTCTAGCAGAGTATTAAT
TACCACTGACCTGTTGGCCAGAGGCATTGATGTGCAGCAGGTCTCCCTAGTCATCAATTATGACCTTCCC
ACCAACAGGGAAAACATCCACAGAATTGGTCGAGGTGGTCCGTTTGGCCGAAAAGGTGTGGCTATTA
ACATGGTGACAGAAGAAGACAAGAGGACTCTTCGAGACATTGAGACTTTCTACAACACCTCCATTGAAGA
GATGCCCTCAACGTTGCTGACCTCATT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	1218 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_199372.2 , NP_955404.1
RefSeq Size:	1801 bp
RefSeq ORF:	1221 bp
Locus ID:	287436
Cytogenetics:	10q24
MW:	46.2 kDa
Gene Summary:	mouse homolog is a protein synthesis initiation factor [RGD, Feb 2006]