

## Product datasheet for **RR211389L4V**

### Furin (NM\_019331) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Furin (NM_019331) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Furin
Synonyms:	Pace; Pcsk3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019331
ORF Size:	2379 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211389).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019331.1</a> , <a href="#">NP_062204.1</a>
RefSeq Size:	4259 bp
RefSeq ORF:	2382 bp
Locus ID:	54281
UniProt ID:	<a href="#">P23377</a>
Cytogenetics:	1q31
Gene Summary:	may be one of several enzymes involved in the proteolytic processing of extracellular superoxide dismutase (EC-SOD) [RGD, Feb 2006]



[View online »](#)