

Product datasheet for **RR211408L3V**

Usp11 (NM_001008861) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp11 (NM_001008861) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Usp11
Synonyms:	Uhx1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008861
ORF Size:	2763 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211408).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008861.2 , NP_001008861.2
RefSeq Size:	3423 bp
RefSeq ORF:	2766 bp
Locus ID:	408217
UniProt ID:	Q5D006
Cytogenetics:	Xq11
Gene Summary:	deubiquitinating enzyme; cleaves ubiquitin from ubiquitin-conjugated protein substrates; important for processes such as cell cycle progression and transcriptional activation [RGD, Feb 2006]



[View online »](#)