

## Product datasheet for **RR211433**

### **Pcbp2 (NM\_001013223) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pcbp2 (NM\_001013223) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Pcbp2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR211433 representing NM\_001013223  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGACACCGGTGTGATTGAAGGTGGATTAATATGTCCTCACCATCAGGCTACTTATGCACGGAAAGG  
AAGTTGGCAGTATCATCGAAAGAAAGGAGAATCTGTTAAGAAGATGCGCGAGGAGAGTGGTGCACGTAT  
CAACATCTCAGAAGGGAATTGCTCTGAGAGAATTATAACTTTGGCTGGACCGACTAATGCCATCTTCAA  
GCCTTTGCTATGATCATTGACAACTGGAAGAGGACATAAGCAGCTCTATGACCAATAGCACAGCTGCCA  
GTAGACCCCGGTTACCTTACGGCTGGTGGTTCCTGCTAGTCAGTGTGGCTCTCTTATTGAAAAGGTGG  
TTGCAAGATCAAGGAAATAAGAGAGAGTACAGGGGCTCAGGTCCAGGTGGCAGGGGATATGCTCCCCAAC  
TCAACTGAGCGAGCAATCACTATTGCTGGCATTCCGCAGTCCATCATTGAGTGTGTCAAACAGATCTGCG  
TGGTCATGTTGGAGACTCTCTCCAGTCCCCCGAAGGGCGTGACCATCCCGTACCGGCCAACGCCCTC  
CAGTTCTCCAGTCATCTTTGCAGGTGGTCAGGACAGGTACAGCACAGGCAGCGACAGTGGCAGCTTTCC  
CACACCACCCGTCATGTGCTCAACCCTGACCTGGAGGGACCCTCTAGAGGCTATACCATTCAAG  
GACAGTATGCCATTCCACAGCCAGATTTGACCAAGCTGCACCAGTTGGCAATGCAACAGTCTCATTTTCC  
CATGACCCATGGCAACACCGGATTCAAGTGGCATTGAATCCAGCTCTCAGAGGTGAAAGGCTATTGGGGT  
TTGGATGCATCTGCTCAAACACTTCTCATGAACTCACCATTCCAAATGATTTGATTGGTGCATATTG  
GGCGTCAGGGCGCCAAAATCAATGAGATCCGTCAGATGTCTGGGGCGCAGATCAAATGCAAAACCCAGT  
GGAAGGATCTACTGATCGGCAGGTTACCATCACTGGATCTGCTGCCAGCATTAGTCTGGCTCAGTATCTA  
ATCAATGTCAGGCTTTCCTCGGAGACAGGTGGCATGGGGAGCAGC

**ACGCGT**ACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR211433 representing NM\_001013223  
 Red=Cloning site Green=Tags(s)

MDTGVIIEGGLNVTLTIRLLMHGKEVGSIIIGKKGESVKKMREESGARINISEGNCPERIITLAGPTNAIFK  
 AFAMIIDKLEEDISSSMTNSTAASRPVTLRLVVPASQCGSLIGKGGCKIKEIRESTGAQVQVAGDMLPN  
 STERAITIAGIPQSIIECVKQICVVMLETLSSQSPPKGVTIPYRPKPSSSPVIFAGQDRYSTGSDSASFP  
 HTTPSMCLNPDLEGPPLEAYTIQGYAIPQDLTKLHQLAMQQSHFPMTHGNTGFSGIESSPEVKGYWG  
 LDASAQTTSHELTIPNDLIGCIIGRQGAKINEIRQMSGAQIKIANPVVEGSTDRQVTITGSAASISLAQYL  
 INVRLSSETGGMGSS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001013223

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORF Size:</b>              | 1095 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001013223.1</a> , <a href="#">NP_001013241.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>           | 1712 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1098 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 363005                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 7q36                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>MW:</b>                    | 38.6 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |