

Product datasheet for **RR211480**

Ncbp3 (NM_001108281) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ncbp3 (NM_001108281) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ncbp3
Synonyms:	RGD1308139
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR211480 representing NM_001108281
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGCGCTTCGCGGCTCGCAGTGTCTGTGAAAGCGAGGCCCGCTGGGCGGCTTTGGGGCTTC
 CGTCGCCCCGAGGTGGAATCCGGTTTGGAGCGTGGCGAGCCGGAACCCATGGAGGTTGAAGAAGCGAGCT
 GGAAATCGTGCCGGTGGCGGCTCCCTAAAGGAACTGCTCCCGGACACAAGCAGAAGATATGAAAAAAG
 GCCGCGAGCTTCACTACTGGGATCGATGTACCTCCAAGGAAGCAATTGAAAAGAAAAGAACGCGAGCCA
 AACGCTTCCACTTTGAGCAGAAGTAACTTGGCCAGAGGAATGTAGCCTTGGACAGAGACATGATGAA
 GAAAGCAATCCCCAAGTAAGACTGGAGACGATCTACATTTGTGGCGTGGATGAGATGAGTACCCAGGAC
 ATCTTTTCTACTTTAAGGAGTACCCTCCAGCTCACATTGAGTGGCTGGATGACACCTCTGTAATGTGG
 TTTGGCTAGATGAGATGACAGCCACAGGCTCTCATCAATATGAGTTCTCTGCCTGCCAGGATAAGAT
 AAGAAGCAGAGAAGCCAGTGAAGACAAGTCGTCTGAGAAAAGACAAGAAAGACAAGCAGGAAGACAGTTCA
 GATGACGACGAGACTGAAGAAGGAGAAGTTGAAGATGAGAATCAAGTGTAGAGTTGGACACATTGT
 CTCAGGTAGAAGAAGAGTCTCTGTTAAGAAATGACCTTCGTCCAGCTAACAACTCGCTAAAGGAAACAG
 ACTGTTTATGAGATTTGCTACGAAAAGTGAACAAAAGGAACTTGGAGCAGCCAGAAGAAGTCAATATTAC
 ATGAAATATGAAATCCAAATATGGAGGCATGAAAGGAATTCGAGCAATTCATGGAAGCGAAGATATC
 ATCCCGTCGCATTACGCGGGATGTGATCAAGAAGAGAGCCCTGATTGGGGATGACGTTGGCTTGACGTC
 GTATAACACCGACATCCGGATTAGTAAACGTTCTGAGGAACCCATTGAGGAAGAGGAGGAGGAAGAG
 CAGGAGGAAGAGGACCAGGACATGGATGCTGATGACAGAGTCGTGGTGGAGTACCATGAGGAGCTCCAG
 CTCTCAAGCAGCCCCGTGAGCGGAGCATGTCTAGGCATCCAGTGGCAGCAGCTCCGACTCCGACGAGAT
 GGACTATGATCTGGAAGTGAAGTATCTCCACTCCCTCCCAAAGAAAAGCATGAAGATGACCATGTAT
 GCAGATGAAGTTGAGTCTCAGTTGAAAAGCATTAGGAATCCCATGAGGTCAGATAGTATATCCACAAGCA
 ATATAAAAAACGAATTGGTAACAAATGCCTCCTGAGAAGTTTGCAGATGTCCGGCACCTTTTAGAGGA
 AAAGCGTCAGCATTCTGTCCACGGCCAGCAGGCGGCAGTACTAAACCAGATATTCGCCAGCGGTTAGGA
 AAAAGACCGTATTCTCCAGAAAAGGCTTTTAGCAGTAACCAGGTGGTCCGGAGAGAGCCGCTTTCGGATG
 TACATAGTAGGCTAGGTGTTCCAGGCAAGATGTTAAAGGCTCTACTCTGATACTCGGGAGAGGAAGTC
 AGGTGGTTTATGGACTCGCTTAGGATCTACACCAAGACCAAAGAAAAGAACACGAAGAAGTGGATCAT
 AGGGCATCTGGGACTGAGGAAGTGAAGTCTCAGAGCTTCCAGGGCATGGGGTCTCTGATTAAAGAGAAA
 AACAGTCTCGACAGAAGAAGAGTCGGCTGGATAGCTTCCCGTCACTCCAGATCGAAGTACGCGGGAAG
 TAGCTCTGGTTCAGAGGCAGAGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RR211480 representing NM_001108281
 Red=Cloning site Green=Tags(s)

MAVRGLRVSVKAEAPAGPALGLPSPEVESGLERGEPEPEMEVEEGELEIVPVRRSLKELLPDTSRRYENK
 AGSFITGIDVTSKEAIEKKEQRAKRFHFRAEVNLAQRNVALDRDMMKKAIPKVRLETIYICGVDEMSTQD
 IFSYFKEYPPAHIEWLDDTSCNVVWLDEMATRALINMSSLPAQDKIRSREASEDKSSEKDKKQEDSS
 DDETEEGEVEDENSSDVELDTLSQVEEESLLRNDRPANKLAKGNRFLMRFATKDDKELGAARRSQYY
 MKYGNPNYGGMGILSNSWKRRYHSRRIQRDVIKKRALIGDDVGLTSYKHRHSGLVNVPPEPIEEEEEE
 QEEEDQMDADDRVVVEYHEELPALKQPRERSMSRRSSASSSDSDMDYDLELRMISTPSPKSMKMTMY
 ADEVESQLKSIRNPMRSDSISTSNIKNRIGNKLPPEKFADVRHLLLEEKQHSRPRPAGGSTKPDIRQLG
 KRYPYSPEKAFSSNQVVRREPSSDVHSRLGVPRQDVKGLYSDTRERKSGGLWTRLGSTPKTKEKNTKKVDH
 RASGTEEDDSELQRWALIKEKEQSRQKSRDLPSLQIEVSRESSSGSEAES


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108281

ORF Size: 1845 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001108281.1 , NP_001101751.1
RefSeq Size:	2880 bp
RefSeq ORF:	1848 bp
Locus ID:	360563
Cytogenetics:	10q24
MW:	70.1 kDa