

Product datasheet for **RR211487**

Mtmr7 (NM_001107312) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mtmr7 (NM_001107312) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mtmr7
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR211487 representing NM_001107312 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAACACATCCGTACCCCTAAGATACTTCACAGTCAGATCTCCACCATTGAAAAGCAGGGCAGCAGACAG
CTACGGGGTCCCTCTACTCATTTCGCTGCAAGAAGCTTTCAGATTGTACAACCTGTACATCCGCAGGAGAG
AGACTGCCACGAGTGTATGTCTCTGATACGCCTTGCAGGGCCAGTGAGATACGAGGAGCTATACTGC
TTTTCTTCAACCCCAAGCTGGATAAGGAAGAACGAGAGCAGGGCTGGTTGCTGGTTGACCTGAGTGAAG
AATACAAGCGAATGGGCCTCCCGGATAATTATTGGCAGCTCAGCGATGTGAATAGAGACTACAGAGTTTG
TGACTCTTACCCACGGAAGTGTATGTCCCGAGATCGGCCACAGCACACATCATCGTGGCAGCTCCAAG
TTCCGGAGCAGACGGCGATTTCTGCCCTTTCCTACTACTGCAAAGAAAGCCATGCCTCCATCTGCAGGA
GCAGCCAGCCCTTTCGGCTTCAGTGTCTGGTGCCTGGAGGATGAGCAGATGCTCCAAGCCATCAGGAA
AGCCAATCCCAGAAGTACTTCATTTACGTCGTGGACACGAGGCCTAACTTAATGCAATGGCAAATCGT
GCTGCAGGAAAGGCTATGAGAATGAAGACAATTATTCCAATATCAAGTTTCAGTTTATTGGGATAGAGA
ACATCCACGTCATGAGGAACAGTCTTCAGAAAATGCTGGAAGTATGTGAACCTAAATCTCCATCCATGAG
TGATTTCTATGGGGCTCGAGAACTCAGGCTGGCTGAGGCATATCAAAGCCATCATGGATGCAGGGATC
TTCAATTGCAAAGGCAGTTTCAGAGGAAGGGACAAGTGTGCTTGTCCACTGTTCTGACGGTTGGGACAGAA
CGGCCAGGTGTCTCGGTGGCAAGCCTGCTGCTGGACCTTACTATCGGACTCTGAAAGGATTCATGGT
GTTAATTGAAAAGGACTGGATTTCTTTGGACATAAAATTAATCATAGATATGGCAATTTAGACGGTGAC
CCAAAAGAAATCTCTCCAGTAATTGACCAGTTCATCGAGTGTGTCTGGCAGTTAATGGAACAGTTCCCTT
GCGCCTTTGAGTTCAATGAGAGGTTTTTGATTACATTCAACATCATGTTTACTCCTGCCAGTTTGGGAA
CTTCTGTGTAATAGCCAGAAGGAGAGACGAGAAGTCAATGAGAATGGAAAAAGGTTTCACAGTTCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR211487 representing NM_001107312
 Red=Cloning site Green=Tags(s)

MEHIRTPKILHSQISTIEKQATTATGCP LLIRCKNFQIVQLVIPQERDCHDVYVSLIRLARPVRYEELYC
 FSNPKLDKEEREQGWLLVDLSEEYKRMGLPDNYWQLSDVNRD YRVCD SYPTEL YVPRSATAHIIVGSSK
 FRSRRRFPALSYCKESHASICRSSQPLSGFSARCLEDEQMLQAIRKANPRSDFIYVVDTRPKLNAMANR
 AAGKGYENEDNYSNIKFQFIGIENIHVMRNSLQKMLEVCELKSPSMSDFLWGLENSGWL RHIKAIMDAGI
 FIAKAVSEEGTSVLVHCS DGWDRTAQVCSVASLLLD PYYRTLKGF MVLEIKDWISFGHKFNHRYGNLDGD
 PKEISPVIDQFIECVWQLMEQFP CAFE FNERFLIHIQH HVYSCQFGNFLCNSQKERRELNENGRVHSS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107312

ORF Size:	1257 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107312.1 , NP_001100782.1
RefSeq Size:	2669 bp
RefSeq ORF:	1260 bp
Locus ID:	306490
Cytogenetics:	16q12.1
MW:	48.7 kDa