

Product datasheet for **RR211567L4V**

Vps4a (NM_145678) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Vps4a (NM_145678) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Vps4a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145678
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211567).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145678.1 , NP_663711.1
RefSeq Size:	1314 bp
RefSeq ORF:	1314 bp
Locus ID:	246772
UniProt ID:	Q793F9
Cytogenetics:	19q12
Gene Summary:	binds Rho family GTPase Rnd2; may play a role in endosomal trafficking [RGD, Feb 2006]



[View online »](#)