

## Product datasheet for RR211568L3

### Smpd3 (NM\_053605) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Smpd3 (NM_053605) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Smpd3
Synonyms:	cca1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211568).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 94338  
UniProt ID: [O35049](#)  
Cytogenetics: 19q12  
Gene Summary: mRNA accumulation occurs in a growth arrested confluent monolayer but not in subconfluent cells; may play a role in regulation of cell growth [RGD, Feb 2006]