

## Product datasheet for **RR211579L4V**

### Prpsap2 (NM\_057131) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Prpsap2 (NM_057131) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Prpsap2
Synonyms:	Pap41
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_057131
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211579).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_057131.1</a> , <a href="#">NP_476472.1</a>
RefSeq Size:	1612 bp
RefSeq ORF:	1110 bp
Locus ID:	117272
UniProt ID:	<a href="#">O08618</a>
Cytogenetics:	10q22
Gene Summary:	associates with the catalytic subunits and 39 kDa associated protein of the ribose-phosphate diphosphokinase enzyme complex [RGD, Feb 2006]



[View online »](#)