

## Product datasheet for **RR211603L3V**

### Wfdc1 (NM\_133581) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Wfdc1 (NM_133581) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Wfdc1
Synonyms:	ps20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133581
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211603).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133581.1</a> , <a href="#">NP_598265.1</a>
RefSeq Size:	1071 bp
RefSeq ORF:	639 bp
Locus ID:	171112
UniProt ID:	<a href="#">O70280</a>
Cytogenetics:	19q12
Gene Summary:	involved in the regulation of growth [RGD, Feb 2006]



[View online »](#)