

## Product datasheet for **RR211604L3V**

### Gap43 (NM\_017195) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gap43 (NM_017195) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gap43
Synonyms:	Basp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017195
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211604).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017195.3</a> , <a href="#">NP_058891.1</a>
RefSeq Size:	1499 bp
RefSeq ORF:	681 bp
Locus ID:	29423
UniProt ID:	<a href="#">P07936</a>
Cytogenetics:	11q21
Gene Summary:	may play a role in neurite outgrowth and nerve regeneration [RGD, Feb 2006]



[View online »](#)