

Product datasheet for **RR211618**

Xrcc2 (NM_001109215) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Xrcc2 (NM_001109215) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Xrcc2
Synonyms: RGD1564823
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR211618 representing NM_001109215
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCAGCGACTTTCGAGGGCAGAATCCGGGACCGAGCTCCTTGCCCGACTTGAAGGCAGAAGCTCCT
TAAAAGAACTAGAGCCCAGCCTGTTTGCTGATGAAGATTCACCAGTGCACGGTGATATTCTTGAATTTCA
TGGTCCAGAAGGAACGGGAAAACAGAAATGCTTTATCATTTAACAGCCGATGTATACTCCAAAATCA
GAGGGTGGACTGCAAATAGAAGTCTTATTTATTGATACAGATTACCACTTTGACATGCTACGGCTTGTGA
CAGTTCTTGAGCACAGACTCTCGCAGAGCTCGGAGGAAGCCATCAAGCTGTGCCTGGGCAGGCTGTTCTCT
GGCGTACTGCAGTAGCAGCATGCAGCTGCTGCTCACGCTGTAAGGCTGTAAGGCTGTAAGGCTGTAAGGCT
CCCTCTCTGCTGCTCATTGTGGATAGCATGTCATCTTTCTACTGGATAGACCGGGTCACTGGAGGAG
AGAGTGTGAGCCTGCAGGAGTCAACACTGCAGAAGTCTCTCAGCTCCTGGAGAGGCTCGTCACCGAGTA
CCGCTTGCTGCTGTTTCGCAACAACACAGAGTCTGTTGCAGAAAGGCTCAGACTCAGCAGAAGGGCTTCC
TCCTCCAAGCACCCGTGCGATGGAGACATGGACTACAGAGCCTATCTCTGCAAGGCTGGCAGAGGGTTG
TGAAGCACAGAGTCGCTTCTCCAGAGAGGATGAGGCCAAGAGCAGCCGGTTCGCACTAGTTTCACGACA
TTAAAAAGTAACAGTTTAAAAAACACGTTTTTATGATTAGAGAAAGTGGGGTGGAGTTTTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR211618 representing NM_001109215
Red=Cloning site Green=Tags(s)

MCSDFRRAESGTELLARLEGRSSLKELEPSLFADEDSPVHGDILEFHGPEGTGKTEMLYHLTARCILPKS
 EGGLEIEVLFIDTDYHFDMLRLVTVLEHRLSQSSEEAIKLCGLFLAYCSSMQLLLTLYSLEALLCSH
 PSLCLLIVDSMSSFYWIDRVSGGESVSLQESTLQKCSQLLERLVTEYRLLLFATTQSLQKGSDSAEGPS
 SSKHPCDGDMDYRAYLCKAWQRVVKHRVVF SREDEAKSSRFALVSRHLKSNLSLKKHVFMIRESGVEFC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001109215
ORF Size: 834 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109215.1 , NP_001102685.1
RefSeq Size:	3509 bp
RefSeq ORF:	837 bp
Locus ID:	499966
Cytogenetics:	4q11
MW:	31.4 kDa