

Product datasheet for **RR211623L4V**

Tmem126b (NM_001106280) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem126b (NM_001106280) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem126b
Synonyms:	RGD1308371
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001106280
ORF Size:	687 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211623).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001106280.2 , NP_001099750.2
RefSeq Size:	870 bp
RefSeq ORF:	690 bp
Locus ID:	293114
UniProt ID:	B2RZD2
Cytogenetics:	1q32
Gene Summary:	Chaperone protein involved in the assembly of the mitochondrial NADH:ubiquinone oxidoreductase complex (complex I). Participates in constructing the membrane arm of complex I (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)