

Product datasheet for RR211641

Meox1 (NM_001108837) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Meox1 (NM_001108837) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Meox1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR211641 representing NM_001108837
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGATCCGGTGGTCAACAGCTGCGTGAGGAACCCAGCCCCAGCTCCCGTCTGGGGCTGCCTTCGAA
 ACCCCCACTCGGAAGGCAGCAGCGCCTCAGGGCTGCCCCATTACCCACCAACCCCGTTTTCTTCCACCA
 AAAGTCCGATTTCCAGCGACAGCAGCATACCCCGACTTCTACTTCTGCTGGCAGCCACCCACAC
 AGCCTGCCCGCGCTGAGCGAATCTTCAATGAGCAGCATCCTGCCTTCCACAGACCCCTGACTGGCATT
 TCCCCATCTCTGAAGCGGGCAAAGGCTCAACCTAGGCCAGCTGGGAGCGCCAGAGAGATGGGGCTGG
 CAGCCCCGGCCTGGTGGACGGCACAGGAGGATTGGGTGAGGACTGCATGGTGTGGGACCATCGCCAT
 GAGACGGAGAAGAAATTATCCAGAAGGAAAAAGGAGAGGTGAGACAACCCAGAGAATGGAGGAGAAAGC
 CAGAAGGCAGCAGCAAGGCCCGGAAGGAAAGGACAGCTTTCACCAAGGAGCAGCTACGGGAGCTAGAGGC
 TGAGTTCGCCACCATAACTACCTGACCAGGCTCCGTAGATACGAGATTGCAGTCAACCTGGACCTTCA
 GAGCGGCAGGTCAAAGTGTGGTCCAGAACCGGAGAATGAAGTGAAGCGTGTGAAGGGGGTCAAGCCG
 TGTCACACAGGAGCAGGACCCAGAGGATGGAGACTCTGCAGCTTCTCCGAGTTCAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR211641 representing NM_001108837
 Red=Cloning site Green=Tags(s)

MDPVVNSCVRNPPAPVWGCLRNPHEGSSASGLPHYPTPF SFHQKSDFPATAAYPDFSTSLAATPH
 SLPRAERIFNEQHPAFPQTPDWHFPISEAGQRLNLGPAGSAREMGAGSPGLVDGTGGLGEDCMVLGTIAH
 ETEKKLSRRKKERSDNPENGGGKPEGSSKARKERTAFTEQLRELEAEFAHNNYLTRLRRYEIAVNLDL
 ERQVKVWFQNRMMKWKRVKGGQPVSPQEQDPEDGDSAASPSSE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108837

ORF Size: 759 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108837.1, NP_001102307.1</u>
RefSeq Size:	1871 bp
RefSeq ORF:	762 bp
Locus ID:	363684
Cytogenetics:	10q32.1
MW:	27.9 kDa