

Product datasheet for **RR211658L4V**

Rpl17 (NM_201415) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpl17 (NM_201415) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rpl17
Synonyms:	AS1; rpL23
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_201415
ORF Size:	552 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211658).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_201415.1 , NP_958818.1
RefSeq Size:	631 bp
RefSeq ORF:	555 bp
Locus ID:	291434
UniProt ID:	P24049
Cytogenetics:	18q12.2
Gene Summary:	component of the 60S subunit of the ribosome, the organelles that catalyze protein synthesis; member of the L22P family of ribosomal proteins; induced under amino acid deprivation of cells [RGD, Feb 2006]



[View online »](#)