

## Product datasheet for **RR211674L4V**

### Lhcgr (NM\_012978) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lhcgr (NM_012978) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lhcgr
Synonyms:	LH/CG-R; Lhr; LSH-R; LSHR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012978
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211674).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012978.1</a> , <a href="#">NP_037110.1</a>
RefSeq Size:	2902 bp
RefSeq ORF:	2103 bp
Locus ID:	25477
UniProt ID:	<a href="#">P16235</a>
Cytogenetics:	6q12
Gene Summary:	for both luteinizing hormone and choriogonadotropin; involved in reproductive development and function [RGD, Feb 2006]



[View online »](#)