

## Product datasheet for **RR211679L3V**

### Pat11 (NM\_001108520) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pat11 (NM_001108520) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pat11
Synonyms:	Pat1b; RGD1305514
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001108520
ORF Size:	2310 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211679).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001108520.2</a> , <a href="#">NP_001101990.1</a>
RefSeq Size:	4007 bp
RefSeq ORF:	2313 bp
Locus ID:	361736
UniProt ID:	<a href="#">B5DF93</a>
Cytogenetics:	1q43



[View online »](#)

**Gene Summary:**

RNA-binding protein involved in deadenylation-dependent decapping of mRNAs, leading to the degradation of mRNAs. Acts as a scaffold protein that connects deadenylation and decapping machinery. Required for cytoplasmic mRNA processing body (P-body) assembly (By similarity).[UniProtKB/Swiss-Prot Function]