

Product datasheet for **RR211690**

Clic6 (NM_176078) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Clic6 (NM_176078) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Clic6
Synonyms:	Clic6b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR211690 representing NM_176078
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCTGAGGCCACGGAGCCCAAGGAGGTCTTTCGGGTTCCCAAGGGCAACCAGAAGGTGCGGTGATCG
 AGGGACCAGGAGAGCCTGGGGCTGCGGACCTGGAGGGCGGGAAGCGAGCGAGGGAGCTGCAGAGGCACC
 CAGGGACCTGGGGGAGGCGCGGAGGCCATGGCCTCAGGGAAGGAGGAAGGAGGCTGTGGCAGGATGGG
 GAAATCGGGGAAGTGCAGGCGCAGGACCCGCGCCCGGGCCCGGGACGGAGACACCGGGCACCAGTGGAG
 CTCCGGGTGAAGCGGAGGCTGCGGAATGCGACTCCGAGGGCGCCTTGATCCCGCAGAGCGGAGGAGGCGC
 CAAGCGCCAGCAGGTGCAGGGAACGAGCTCGGGGCTAGACGCTCAGGGTGAGGCCCTGAGGTCCCCGAG
 GACGCTCGCCGGGAGCCAGAGGACCCCAAGGCTTCAGAGGCAGGGGAAGAGGCAGAATCTGGTCAGGAAG
 CTCTGGGAGGCAGCGGCCAGAGTCGCAGATCAATCCAGAAGTCCAGGGACCTGTAGGGGACAACATGGA
 CACAGAGGCCCCCGCAGGAGAGCCACAGGGGTAGAGGGAGAGCCCCAGGGTGGAGGGGAAAGCAGCCCT
 CAGCCCCAGGATGAGGCGATTGAGATTGACGCCGAGAGGTTGGAGGGCACGAGCCCGGGGAGCTGGCGG
 GCGCTTCAGCCGAGACGCCAAGGGGGAAGGGGAAACCCTGAGGAAAGATGGGTTTCGAGGAAGCTGCTCC
 CGAGGAGGCAAGAGTAGACTCTGGCGAGAATGGGTTTCGAGGAAGCAGCTCCCGAGGAGGCAAGAGTAGAC
 TCTGGAGAGAATAGGGACCAGGGCAGGCTCCAGGAGGAGACCGGGGAGGAAGAAGCGAGACCGGAATCCG
 GGTTGAAGGGACCCTGCGAGGAGGCGATCCAAGAGAAGGCCCCGACGGCAGCCTTGATGGAGAGGAGGC
 CAAATCCACCGGGCAGGAGTGCAGAGTGGAGCTCAGTAACCACTAGCGGAGGAGACCAGCGCCAG
 GCGGAGAAGAGTTGGGGCAGTGAATGGACGTAGGGAGAAGCGCCAGCTTCGAGGAGGGGGATCTGG
 GACAGGAGCACGACATCACCTCTTTGTCAAGGCTGGCTCTGATGGAGAGAGCATCGGAACTGCCATT
 CTCACAGCGTCTCTTTATGATCCTCTGGCTAAAGGGTGTATATTCAATGTGACAACGGTGGACCTGAAA
 CGGAAGCCTGCAGACCTGCAGAACCTGGCTCCTGGGACCAACCCTCCTTTCATGACTTTTGATGTTGAAG
 TCAAGACAGATGTGAATAAGATTGAGGAGTCTTAGAGGAGAAGTTAGTCCCCCGAGGTACCCTAAGCT
 GGGAACCCAACACCTGAATCCAATCAGCGGGAATGACGTTTTTGCCAAATTCTCAGCGTTTATTAAG
 AACACCAAGAAGGATGCAATGACATTTATGAAAAGAACCTACTACGGGCCCTGAAGAAGCTGGACAGTT
 ACTTAACAGCCCGCTGCCAGACGAAATAGATGCCTACGTACGGAGGATGTACTGTTTCCAGAGGAA
 GTTCTCGATGGGGATGAGCTCACACTGGCCGACTGTAACCTCTTACCAAGCTCCACATCATTAAAGATT
 GTGGCCAAGAAATACAGAGGTTTGAAGTTCTTCTGAAATGACTGGCATTGCGGCTACTTGAACAATG
 CTTATGCCCGGACGAGTTCACTAACACGTGTCCAGCTGACCAGGAGATTGAACACGCGTATTACAGACG
 AGCGAAAAGGATGAAG

ACCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RR211690 representing NM_176078
 Red=Cloning site Green=Tags(s)

MAEATEPKVSSGSQQQPEGAVIEGPGEPGAADLEGREASEGAAEAPRDLGEGAEAMASGKEEGCGQDG
 EIGEVQAQDPRPGTETPGTSGAPGEAEAAECDSEGALIPQSGGAKRQVQGTSSGLDAQGEAPEVPE
 DARREPEDPKASEAGEEASGQALGGSAPESQINPEVQGPVGDNMDTEAPAGEPQGSEGEPPQGGESSP
 QPQDEAIEIAAAEVGGHEPGELAGASAADAKGEGETLRKDGFEAAPEEARVDSGENGFEEAAPEEARVD
 SGENRDQGRLEETGEEEARPESGLKGPCEEAIQEKAPDGSLDGEEAKSTGHEESQVELSNHLAEETSAQ
 GGEELGRVNGRRENGPASEEGDLGQEHITLQVKAAGSDGESIGNCPFSQRLFMILWLKGVIFNVTTVDLK
 RKPADLQNLAPGTNPPFMTFDGEVKTVDNKKIEEFLEEKLVPPRYPKLGTQHPESNSAGNDVFAKFSAFIK
 NTKKDANDIYEKNLLRALKKLDSYLNPLPDEIDAYSTEDVTVSQRKFLDDELTLADCNLLPKLHIKI
 VAKKYRGFEFPSEMTGIWRYLNNAYARDEFTNTCPADQEIHAHSDAAKRMK

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_176078

ORF Size: 1836 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_176078.2](#), [NP_788267.1](#)

RefSeq Size: 2252 bp

RefSeq ORF: 1839 bp

Locus ID: 304081


UniProt ID: [Q811Q2](#)

Cytogenetics: 11q11

MW: 64.7 kDa

Gene Summary: putative intracellular chloride channel; binds to dopamine receptors [RGD, Feb 2006]

Product images:



Circular map for RR211690