

Product datasheet for **RR211696**

Tectb (NM_001106200) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tectb (NM_001106200) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tectb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR211696 representing NM_001106200
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGGCCAGGGCCTTCGTTCTGCTGGCCCTCTTTGCGGAAGCTTCGGCGAAATCATGCACTCCGAATA
AAGCAGATGTCATCCTTGTGTTTTGTTACCCCAAGACCATCATCACTAAAATCCCCGAGTGTCCCTACGG
ATGGGAAGTACACCAGCTGGCCCTCGGGGACTGTGTTACAATGGGGTCCATGAAGTGGCTACTACCAA
TTCGTCATCCCAGATCTGTCACCCAAGAACAAGTCCTACTGTGGAACCCAGTCAGAGTACAAGCCTCCCA
TCTACCACTTCTATAGCCACATCGTGTCCAATGACAGCACCGTGATTGTGAAGAACCAGCCTGTCAACTA
CTCCTTCTCCTGCACCTACCACTCCACCTACTTGGTGAACCAGGCTGCCTTTGACCAGAGAGTGGCCACT
GTTTCATGTCAAGAACGGGAGCATGGGCACATTTGAAAGCCAGTTGTCCCTCAACTTCTACACTAATGCCA
AGTTTTCCACCAAGAAAGCTCCCTTCGTCCTGAAACATCTGAAATCGGTTCCGGATCTGTTTGCAGG
AGTAGAAGCCAAGGGACTAAGTGTTCGGTTTAAAGTAGTCTTGAATAGCTGTTGGCCACCCCTCGGCT
GACTTCATGTACCCCTTGCAATGGCAGCTGATCAATAAGGGTTGCTCCACTGATGAAACCGTCTCGTGC
ACGAGAATGGCAAGGACCACAGGGCCACTTCCAATTCAATGCCTTCCGGTTCAGAACATCCCCAACT
ATCCAAGTTTGGTTACTGTGAGACGTTTCATCTGCGACAGTGAGAAGCTCTCCTGCCAGTGAAGTGT
GACAAAACGGAAACGCATGCTACGTGACCAGACAGGGGGTGTCTGTTGTGGAGCTGTCCCTACGGAGCA
GAGGGTTTTCCGGTCTCTGTGACTTCTCAGATGTTCTTCTCCACCTCGTCTGATGCTGGGGACCTGGGC
CGTGTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR211696 representing NM_001106200
 Red=Cloning site Green=Tags(s)

MVARAFVLLALFAEASAKSCTPNKADVILVFCYPKTIITKIPECPYGWEVHQLALGGLCYNGVHEGGYYQ
 FVIPDLSPKNKSYCGTQSEYKPPYHFYSHIVSNDSTVI VKNQPVNYSFCTYHSTYLVNQAAFQDQRVAT
 VHVKNKSMGTFESQLSLNFYTNAKFSTKKEAPFVLETSEIGSDLFAGVEAKGLSVRFKVVNLNSCWATPSA
 DFMYPQLWQLINKGCSTDETVLVHENGKDHRAFQFNAFRFQNI PKLSKVWLHCETFICDSEKLSCPVNC
 DKRKRMLRDQTTGGVLLVLESLRSRGSGLCDFSDVLLHLVLMGTWAVL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106200

ORF Size:	987 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106200.1 , NP_001099670.1
RefSeq Size:	2749 bp
RefSeq ORF:	990 bp
Locus ID:	292124
Cytogenetics:	1q55
MW:	36.9 kDa