

## Product datasheet for **RR211705**

### Scarf2 (NM\_001105864) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Scarf2 (NM_001105864) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Scarf2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**


>RR211705 representing NM\_001105864  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGAGGGCGCAGGGTCCCGGGGGCGGGCCGGCGCGGCCAGGGAGCCGAGGGCCCGCTGCTGC  
 TGCTGCTCTGGCTGCTGCCCGTCTCGCGCGGCCAGGACTTGAACCCGCGAGGCCGAACGTGTGCCG  
 CACGCCCGTTCTCAGGTGCTCAGTGTGCTGCTGGTGGAGGCAGCTGGGGGATGAGTGTGGGATCGCG  
 GTGTGTGAAGCAATCCACGTGCTCAGAAAATGAGGTGTGCGTGCGGCCAGGCGAGTGCCGCTGCCGAC  
 ATGGCTACTTCGGTGCCAACTGCGACACTAAGTGCCCGCGCCAGTTCTGGGGCCCGACTGCAAGGAGAG  
 GTGCAGTTGCCATCCGCACGGACAGTGCAGGACGTGACCGGACAGTGTACGTGTACGCGCGCCGCTGG  
 GGTGCGCGCTGCGAGCATGCTTGCCAGTGCCAGCATGGCAGTGCACCCCGGGAGCGCGCGTGGCCGT  
 GCGAGCCGGGCTGGTGGGTGCACAATGCGCTAGCGCCTGCTACTGCAGCGCTACTTCGCGCTGCGATCC  
 ACAGACAGGCGCCTGCCTGTGCCACGCAGGCTGGTGGGGCCGAGCTGCAACAACAGTGTCTCTGCAAC  
 TCGTCGCCCTGCGAGCAGCAGAGCGGCCGCTGCCAGTGTGCGGAGCGCAGTTTCGGCGCGCGCTGTGATC  
 GCTATTGCCAGTCTCCGCGGTGCTGCCACCCGGTGGACGGCACGTGCGCCTGCGACCCGGGCTACCG  
 GGGCAAGTATTGTCGTGAGCCATGCCCGCTGGCTTCTACGGCCTCGGGTGTGCGCGTGGTGTGGCCAA  
 TGCAAAGGCCAGCAGCCGTGCACAGTAGTGAAGGCCGCTGTCTGACGTGCAACCCGGCTGGAACGGAA  
 CCAATGTGACCAACCTGCGCCACTGGTTTTATGGCGAGGGTGTGGCCACCCTGTCCACCCTGCCG  
 CGATGGGCATGCCTGCAACCATGTACCCGGCAAGTGTACGCACTGCAATGCGGGTGGATCGCGGACCGG  
 TGCGAGACTAAGTGCAGTAATGGCACTTACGGTGAAGACTGTGCCTTCGTGTGTTCTGACTGCGGCAAGC  
 GCCAGTGTGACTCCAGTCCGGCGCTGCCTTTGCGAGCCCTGGTGTTCACGGACCTACTGTAATGTAAC  
 TTGCCCGGCCGGCTCCATGGTGTGACTGCGCCAGGCCCTGCAGCTGTACAGAGGAATCATGCGACCCG  
 GTCACCCGGGCTGCCACCTGGAGACCAATCAGCGCAAAGGTGTGATGGGCGCCGCGCGCTGCTACGC  
 TGCTCCTCGCGCTGCTGCTCCTGCTGCTCGTTGTTGCTGTGCCTGCCGGGGCAAGGACTCGGCGCGCCG  
 GGAGCTCACCTCGAAGGAAGAAGGCACCTCAACGCTTTTGTGGAAGCTTACGCCGATCAGCATGAAG  
 CTACCCGAATCCCGCTTCGAGGCAGAAGCTGCCCAAAGTCGTAGTGGCCATCATGACTTAGACAACA  
 CGCTCACTGCAGTTTCTGGAACCACCTTCAGGGCTGGAGCAGCCCTCCATCATGGTCTCCAGGGC  
 CTCCTTCTCATATTTGACACCACAGATGAAGGCCCTGTGTAAGTGTGCCCCACGAGGAGGCGACTGCT  
 GAGAACCAGACCCAGAAGCTTCTGCTCTCTCACTGAGGTGCCCGCTGTGTCCCTAGCGCCATGGGCA  
 CTTCCGGTCCCGGGGATGAGGCAGCAGCCCTCCCTGCATCCTCAGACAGCGAGCGCTCAGCATCAAGCGT  
 GGAGGGGCGGGCGGGCGCTGTACGCGCGTGTGGCCCGCGCAGGCCCCGGCCGGCCCGGACCCGATGT  
 GAGGCTGGGGCTTGTCACTGTCCCATCGCCGAGCGCAGGAAGCCACCGCCACCCGACCTGCCACCA  
 AGCCTAAGGTGTCTGGATCCACGGCAAGCAAGGCGCCGTGGCCGCTGCACCTTCGCCCGCCCGCAGG  
 CCGCAAGGCTGCGCCGAGCCCTAGCGGGAGGAAACGGACCCCAAGCAACATGTGGTGCAGCCCCGGGC  
 CTGACCGAAGAAGCCACCGGCCCGGCTCACCCACGCCACCCCGCGCCGAGCGCGCGCCGGCCCGG  
 GTCTCTTGAGCCACGACGCGGGAGGTCCCCGCGCAGTGCGCCGAGGCCCTCTATGCTGGCAGC  
 CGAGCTGCGCGACAAGACTCGCAGCCTGGTTCGCGCCGAGAAGCCCGCCGCTCAGAAGGCCAAGCGC  
 TCCGTGCTCCCGCGTGCACCGCCCGCGGCTTCGGCGTCCGAAGCCTCGGCGTCTGAGAAGGTGGCAG  
 CCAGCGCACCTGAGCCGAGACCTCGGAAGAAGACCCCATCCAGAAACCGCCTCGAAAGAAGAGCAG  
 GGAGGCTGCGGGGAGCCGAGCAGGGCGGGCGGCCCGGAGCTTCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



**Plasmid Map:**


**ACCN:** NM\_001105864

**ORF Size:** 2499 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001105864.1](#), [NP\\_001099334.1](#)

**RefSeq Size:** 3312 bp

**RefSeq ORF:** 2502 bp

**Locus ID:** 287949

**Cytogenetics:** 11q23

**MW:** 88 kDa