

Product datasheet for RR211738L4

Renbp (NM_031095) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Renbp (NM_031095) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Renbp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211738).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

RefSeq ORF: 1260 bp
Locus ID: 81759
UniProt ID: [P51607](#)
Cytogenetics: Xq37
Gene Summary: has N-acetyl-D-glucosamine 2-epimerase activity; inhibits renin activity [RGD, Feb 2006]