

Product datasheet for **RR211767L3V**

Rps6ka1 (NM_031107) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rps6ka1 (NM_031107) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rps6ka1
Synonyms:	S6K-alpha 1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031107
ORF Size:	2205 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211767).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031107.1 , NP_112369.1
RefSeq Size:	3131 bp
RefSeq ORF:	2208 bp
Locus ID:	81771
UniProt ID:	Q63531
Cytogenetics:	5q36
Gene Summary:	catalyzes the phosphorylation of ETS transcription factor ER81 and other transcription factors [RGD, Feb 2006]



[View online »](#)