

## Product datasheet for **RR211770**

### **Bcl2l10 (NM\_053733) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bcl2l10 (NM\_053733) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Bcl2l10  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR211770 representing NM\_053733  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGTGACCCGCTGCAGGATCGCACTAGACGGCTGCTGACTGACTACATATTGTTCTGCGCACGGGCGC  
CGAACACCCTGAGCCACTGCCACGTCTGTTGAGGCGGCCTTGCTGCGCTCTGTGACTAGTCAGATCCA  
ACAGGAGCACCAGGATCTTTTCAACTCCTCCGCGACTACCAGGGCAACCGCTGGAGCTGGTGACACAG  
ATGGCGGATGAGTTGCTCTCCAATGACCAAGAGTTCAACTGGGGCCGCCTGGTGATGCTCCTGGCCTTCG  
TGGGGACGCTAATGAACCAAGACAGGACTGTTAAGCGGAGGAGGATCAAAGAAACCGTCTCCTACTGGA  
GCGAGACTGCTATCTCATAGTGAGCTTGTGTACAATCGACTCACAGGACGGCATCGCTCCTGGCTGGAG  
GCTCACGGTGGCTGGATGGCTTTTGCCAATTCTTCAAGAACCCTTACCACCGGCTTCTGGAGAAGAT  
TGCTGATCCGGGCTATTCTGTCTGTTTCTTGAACGGCCATCTTTTATATCTGGAAATGTTTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR211770 representing NM\_053733  
Red=Cloning site Green=Tags(s)

MGDPLQDRTRRLLDYILFCARAPNTPEPLPTSVEAALLRSVTSQIQEHQDLFNSFRDYQGNRLELVTO  
MADELLSNDQEFNWGRLVMLLAFVGTLMNQDRTVKRRRDQRNRLLLERDCYLIVSLLYNRLTGRHRSWLE  
AHGGWDGFCQFFKNPLPPGFWRRLILIRAILSCFFATAIFYIWKCL

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_053733

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_053733.1</a></u> , <u><a href="#">NP_446185.1</a></u>
<b>RefSeq Size:</b>	1074 bp
<b>RefSeq ORF:</b>	558 bp
<b>Locus ID:</b>	114552
<b>UniProt ID:</b>	<u><a href="#">Q99M66</a></u>
<b>Cytogenetics:</b>	8q24
<b>MW:</b>	21.9 kDa
<b>Gene Summary:</b>	Promotes cell survival. Suppresses apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]