

Product datasheet for **RR211829L4V**

Sidt2 (NM_001108142) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sidt2 (NM_001108142) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sidt2
Synonyms:	RGD1308311
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001108142
ORF Size:	2937 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211829).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001108142.1 , NP_001101612.1
RefSeq Size:	4741 bp
RefSeq ORF:	2940 bp
Locus ID:	315617
UniProt ID:	D3ZEH5
Cytogenetics:	8q22



[View online »](#)

Gene Summary:

Mediates the translocation of RNA and DNA across the lysosomal membrane during RNA and DNA autophagy (RDA), a process in which RNA or DNA is directly imported into lysosomes in an ATP-dependent manner, and degraded. Involved in the uptake of single-stranded oligonucleotides by living cells, a process called gymnosis (By similarity). In vitro, mediates the uptake of linear DNA more efficiently than that of circular DNA, but exhibits similar uptake efficacy toward RNA and DNA. Binds long double-stranded RNA (dsRNA) (500 - 700 base pairs), but not dsRNA shorter than 100 bp (By similarity).[UniProtKB/Swiss-Prot Function]