

## Product datasheet for **RR211841**

### Gspt2 (NM\_001109319) Rat Tagged ORF Clone

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | Expression Plasmids                       |
| Product Name:             | Gspt2 (NM_001109319) Rat Tagged ORF Clone |
| Tag:                      | Myc-DDK                                   |
| Symbol:                   | Gspt2                                     |
| Synonyms:                 | RGD1563213                                |
| Mammalian Cell Selection: | Neomycin                                  |
| Vector:                   | pCMV6-Entry (PS100001)                    |
| E. coli Selection:        | Kanamycin (25 ug/mL)                      |



[View online »](#)

ORF Nucleotide  
Sequence:

>RR211841 representing NM\_001109319  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGATCTCGGCAGCAGCAGCGACTCGGCTCCCGACTGCTGGGATCAGGTGGACATGGAAGCCCCGGTT  
CGGCTCCGAGCGGGATGAAATCGCCCCGCAGCCACGGCAGCGGCTGAGGCCGGGAAGCCGAGGCCCA  
GCGCAGGCATCTCAGCTTGGCCTTCAGCAGTCAGCTCAACATCAACGCCAAACCTTTCGTGCCTAGCGTA  
AGCGCAGCGGAGTTCGTGCCGTCTTCTGCCGGATCGGCCAGCTGCCCGCCCCACAGCCTCCAGCT  
GCGACGAAACCTGCATCGGTGGCGCCGGAGCCTCAAGGTAAAAGGATGGAATGGGAGACCTGTGGA  
ACCTTCCAAGATGGACCTTTAGTGTCTGGGAGGTTCCAATTCAGCTGTTACCATGGAACCTTTCAGAA  
CCTGTTGTAGAAAATGGAGAGGTGGAGATGGCCCTGGAAGAATCGTGGGAGCTTAAAGAAGTAAGTGAAG  
CAGCGCCTGAGGGTTCTTTGGGAGATGTAGGGCCCCAGAAGAGAGTGGCAAGGAAGTATGGAGGAGAA  
AGAGGAAATAAGGAAATCAAATCTGTGGTCAGACCATCAGGTGCTCCTAAGAAAGAACACGTAAATGTG  
GTGTTTCATTGGGCATGTAGATGCCGGCAAGTCAACCATTGGAGGACAGATAATGTTTTTGACAGGAATGG  
TTGACAGAAGGACACTTGAGAAATACGAACGAGAAGCTAAGGAGAAGAATAGAGAGACCTGGTACTTGTG  
CTGGGCCTTAGATACAAACCAGGAAGAAGCAGACAAGGGTAAAACAGTGGAAAGTGGGCCGTGCATATTTT  
GAAACAGAAAAGAAACATTTTCACAATCTTAGATGCCCTGGCCACAAGAGTTTTGTCCCAATATGATTG  
GTGGTGCCTTCTCAAGCTGATTTAGCCGTGCTGGTCATCTCTGCCAGGAAAGGAGAATTTGAAACTGGGT  
TGAAAAAGGTGGCAGACAAGAGAGCATGCAATGTTGGCCAAAACAGCAGGGTGAATACTTAATAGTG  
CTGATTAATAAGATGGATGATCCACAGTAGATTGGAGCAGTGAGCGATACGAAGAATGTAAGAGAAAAC  
TGGTGGCCTTTTTGAAAAAAGTTGGATTCACTCCAAAAAAGGACATTCATTTATGCCCTGCTCAGGACT  
GACTGGAGCAAATATTAAGAGCAATCCGATTTCTGTCCCTGGTACACTGGGTTACCATTTATTCCATAT  
TTGGATAGTTTGCCAACTTCAACAGATCAATTGATGGGCCAATTAGGTTGCCATTGTGGATAAGTACA  
AGGATATGGGCACTGTGGTCTGGGAAAGCTGGAATCAGGGTCCATTTTTAAAGGCCAGCAGCTTGAAT  
GATGCCGAACAAGCACAGTGTGGAAGTTCTAGGAATAGTTCCGATGATGCTGAGACTGATTTTGTAGCT  
CCAGGTGAAAACCTCAAATCAGACTGAAAGGAATCGAAGAAGAAGAGATTCTCCGGGCTTCATACTTT  
GTGAACCAATAATCTCTGCCATTCTGGACGCACATTTGATGTTTCAGATAGTGATCATTGAGCACAAGTC  
TATCATCTGCCAGGTTACAATGCGGTGCTTACATTCATACGTGTATAGAGGAGGTCGAGATAACAGCC  
TTGATCTCCTTGGTAGACAAAAGTCAGGAGAAAAAGCAAGACACGGCCCCGCTTTGTGAAGCAAGATC  
AAGTGTGCATTGCCCGTTAAGGACAGCAGGAATATCTGCCTGGAGACGTTTAAAGATTTTCTCAGAT  
GGTCTGTTTTACTTTAAGAGATGAGGAAAGACCATCGCGATTGAAAAGTTCTGAAGCTGGTCCCAGAG  
AAGGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR211841 representing NM\_001109319  
Red=Cloning site Green=Tags(s)


MDLGSSSDSAPDCWDQVDMEAPGSAPSGDEIAPAATAAAEAAEAQRRLSLAFSSQLNINAKPFVPSV  
 SAAEFVPSFLPGSAQLPAPTASSCDETCIGGAAEPQGRMEWGAPVPEPSKDGPLVSWEGSNSAVTMELSE  
 PVVENGEVEMALEESWELKEVSEAAPEGLGDVGPPEESGKEVMEEKEEIRKSKSVVRPSGAPKKEHVNV  
 VFI GHVDAGKSTIGGQIMFLTGMVDRRTLEKYEREAKEKNRETWYLSWALDTNQEERDKGKTVEVGRAYF  
 ETEKKHFTILDAPGHKSFVPMIGGASQADLAVLVI SARKGEFETGFEKGGQTRHAMLAKTAGVKYLIV  
 LINKMDDPTVDWSSERYEECKELVPFLKKGVSFSPKDIHFMPCSGLTGANIKEQSDFCPWYTG LPIFY  
 LDLSLPNFNRSIDGPIRLPIVDKYKDMGTVVLGKLESGSIFKGGQLVMMPNKSVEVLGIVSDDAETDFVA  
 PGENLKI RLKGI EEEEEILPGFILCEPNNLCHSGRTFDVQIVIIIEHKSIIICPGYNV LHIHTCIEEVEITA  
 LISLVDDKKSGEKSKTRPRFVKQDQVCIARLRTAGTICLETFKDFPQMGRFTLRDEGKTI AIGKVLKLVPE  
 KD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**ACCN:** NM\_001109319

**ORF Size:** 1896 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001109319.1](#), [NP\\_001102789.1](#)

**RefSeq Size:** 2522 bp


**RefSeq ORF:** 1899 bp

**Locus ID:** 501582

**Cytogenetics:** Xq22

**MW:** 69.2 kDa

**Product images:**



Circular map for RR211841