

## Product datasheet for **RR211896L4V**

### Coq6 (NM\_001011983) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Coq6 (NM_001011983) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Coq6
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
ACCN:	NM_001011983
ORF Size:	1428 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211896).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001011983.1</a> , <a href="#">NP_001011983.1</a>
RefSeq Size:	1762 bp
RefSeq ORF:	1431 bp
Locus ID:	299195
Gene Summary:	FAD-dependent monooxygenase required for the C5-ring hydroxylation during ubiquinone biosynthesis. Catalyzes the hydroxylation of 3-polyprenyl-4-hydroxybenzoic acid to 3-polyprenyl-4,5-dihydroxybenzoic acid. The electrons required for the hydroxylation reaction may be funneled indirectly from NADPH via a ferredoxin/ferredoxin reductase system to COQ6.[UniProtKB/Swiss-Prot Function]



[View online »](#)