

Product datasheet for **RR211900L4V**

Stk25 (NM_184049) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stk25 (NM_184049) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Stk25
Synonyms:	MGC94619
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_184049
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211900).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_184049.1 , NP_908938.1
RefSeq Size:	2021 bp
RefSeq ORF:	1281 bp
Locus ID:	373542
Cytogenetics:	9q36
Gene Summary:	human homolog is a member of the Ste20-like kinase family which is activated by cellular, particularly oxidant, stress but which does not activate either of the stress-activated MAP kinase cascades (p38 and SAPKs) [RGD, Feb 2006]



[View online »](#)