

Product datasheet for **RR211903L4V**

Serp1 (NM_030835) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serp1 (NM_030835) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Serp1
Synonyms:	RAMP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030835
ORF Size:	198 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211903).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030835.2 , NP_110462.1
RefSeq Size:	2442 bp
RefSeq ORF:	201 bp
Locus ID:	80881
UniProt ID:	Q9R2C1
Cytogenetics:	2q26
Gene Summary:	may be involved in glycosylation modification of secretory and membrane proteins including glycosylation of MHC class II-associated invariant chain (Ii) [RGD, Feb 2006]



[View online »](#)