

Product datasheet for RR211917L4

Capn3 (NM_017117) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Capn3 (NM_017117) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Capn3
Synonyms:	Lp82; Lp84; Lp85; p94
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211917).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_017117
ORF Size:	2463 bp



[View online »](#)

