

Product datasheet for **RR211918L4V**

Tssc1 (Eipr1) (NM_001012192) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Tssc1 (Eipr1) (NM_001012192) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Eipr1 |
| Synonyms: | Eipr1; Tssc1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001012192 |
| ORF Size: | 1158 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR211918). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001012192.1 , NP_001012192.1 |
| RefSeq Size: | 1603 bp |
| RefSeq ORF: | 1161 bp |
| Locus ID: | 362721 |
| UniProt ID: | Q5PPK9 |
| Cytogenetics: | 6q16 |



[View online »](#)

Gene Summary:

Acts as a component of endosomal retrieval machinery that is involved in protein transport from early endosomes to either recycling endosomes or the trans-Golgi network. Mediates the recruitment of Golgi-associated retrograde protein (GARP) complex to the trans-Golgi network and controls early endosome-to-Golgi transport of internalized protein. Promotes the recycling of internalized transferrin receptor (TFRC) to the plasma membrane through interaction with endosome-associated recycling protein (EARP) complex.[UniProtKB/Swiss-Prot Function]