

## Product datasheet for **RR211934L3V**

### Plk3 (NM\_022187) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Plk3 (NM_022187) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Plk3
Synonyms:	Cnk; Fnk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022187
ORF Size:	1941 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211934).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022187.1</a> , <a href="#">NP_071523.1</a>
RefSeq Size:	2389 bp
RefSeq ORF:	1944 bp
Locus ID:	58936
UniProt ID:	<a href="#">Q9R011</a>
Cytogenetics:	5q36
Gene Summary:	cell cycle associated kinase; may be involved in long-term synaptic plasticity [RGD, Feb 2006]



[View online »](#)