

Product datasheet for **RR211951**

Sox14 (NM_001106850) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sox14 (NM_001106850) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sox14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR211951 representing NM_001106850
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCAAACCTTCAGACCACATCAAGAGGCCCATGAATGCCTTCATGGTATGGTCCCGGGGCCAGCGGC
GCAAGATGGCCCAGGAAAATCCAAGATGCACAACCTCGGAGATCAGCAAACGCCTAGGCGCCGAATGGAA
GCTCCTGTCCGAGGCAGAGAAGCGCCCTACATCGACGAAGCCAAAAGGCTGCGCGCCAGCACATGAAG
GAGCACCCCGACTACAAGTACCGCCCTCGGCGCAAGCCGAAGAACCCTGCTCAAGAAAGACAGGTATGTCT
TCCCCCTGCCCTACCTGGGCGACACGGACCCGCTCAAGGCGGCCGGCCTGCCCGTGGGGGCTCCGACGG
CCTCCTGAGCGGCCCGAGAAAAGCCCGGCCCTTCTTGCCGCGGCCTCGGCGCCTTACTCCTGCTGGAC
CCTGCGCAATTTAGTTCCAGTGCCATCCAGAAGATGGGTGAAGTGCCCCACACGTTGGCCACCAGCGCGC
TGCCCTACGCGTCCACCCTGGGCTATCAGAACGGCGCCTTCGGCAGTCTCAGCTGCCCCAGCCAGCATA
GCACACGCACCCGTCCCCACCAACCCGGGCTACGTGGTGCCCTGTAAGTGTACCGCCTGGTCTGCCTCC
ACCCTGCAGCCCCCTGTCGCTACATCCTTCCAGGCATGACCAAGACTGGCATAGACCCTTATTCGT
CAGCCCACGCTACGGCCATG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR211951 representing NM_001106850
 Red=Cloning site Green=Tags(s)

MSKPSDHIKRPMAFMVWSRQRRKMAQENPKMHNSEISKRLGAEWKLLSEAEKRPYIDEAKRLRAQHMK
 EHPDYKYRPRRKPNLLKKDRYVFPLPYLGTDPLKAAGLPVGASDGLLSAPEKARAFLLPPASAPYSLLD
 PAQFSSSAIQMGVEPHTLATSALPYASTLGYQNGAFGSLSCPSQHTHTHPSPNPGYVVPNCNTAWSAS
 TLQPPVAYILFPGMTKTGIDPYSSAHATAM

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001106850

ORF Size: 720 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106850.1](#), [NP_001100320.1](#)

RefSeq Size: 1937 bp


RefSeq ORF: 723 bp

Locus ID: 300954

Cytogenetics: 8q31

MW: 26.5 kDa

Product images:



Circular map for RR211951