

## Product datasheet for **RR211980L4V**

### Trh (NM\_013046) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Trh (NM_013046) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trh
Synonyms:	Pro-TRH; THR; Trf; TRH01
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013046
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211980).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013046.3</a> , <a href="#">NP_037178.1</a>
RefSeq Size:	1409 bp
RefSeq ORF:	768 bp
Locus ID:	25569
UniProt ID:	<a href="#">P01150</a>
Cytogenetics:	4q34



[View online »](#)

**Gene Summary:**

This gene encodes a member of the thyrotropin-releasing hormone family. Cleavage of the encoded proprotein releases mature thyrotropin-releasing hormone, which is a tripeptide hypothalamic regulatory hormone. The rat proprotein contains five thyrotropin-releasing hormone tripeptides. Thyrotropin-releasing hormone is involved in the regulation and release of thyroid-stimulating hormone, as well as prolactin, and may be involved in cardiac function in heart failure. Disruption of this gene is associated with hypothyroidism, elevated thyroid-stimulating hormone levels, and hyperglycemia. [provided by RefSeq, May 2013]