

Product datasheet for **RR211992L4V**

Gstm4 (NM_001024304) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Gstm4 (NM_001024304) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Gstm4 |
| Synonyms: | MGC108896 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001024304 |
| ORF Size: | 654 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR211992). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001024304.1 , NP_001019475.1 |
| RefSeq Size: | 1583 bp |
| RefSeq ORF: | 657 bp |
| Locus ID: | 499689 |
| UniProt ID: | Q5BK56 |
| Cytogenetics: | 2q34 |



[View online »](#)

Gene Summary:

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Catalyzes the conjugation of leukotriene A4 with reduced glutathione (GSH) to form leukotriene C4. Can also catalyzes the transfer of a glutathionyl group from glutathione (GSH) to 13(S),14(S)-epoxy-docosahexaenoic acid to form maresin conjugate in tissue regeneration 1 (MCTR1), a bioactive lipid mediator that possess potent anti-inflammatory and proresolving actions.[UniProtKB/Swiss-Prot Function]