

Product datasheet for **RR211993L4V**

Insig2 (NM_178091) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Insig2 (NM_178091) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Insig2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178091
ORF Size:	675 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR211993).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178091.4 , NP_835192.1
RefSeq Size:	720 bp
RefSeq ORF:	678 bp
Locus ID:	288985
UniProt ID:	Q80UA9
Cytogenetics:	13q11
Gene Summary:	endoplasmic reticulum protein that blocks proteolytic activation of sterol regulatory element-binding proteins [RGD, Feb 2006]



[View online »](#)