

Product datasheet for **RR212035L4V**

Egln2 (NM_001004083) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Egln2 (NM_001004083) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Egln2
Synonyms:	HIF-PH1; HPH-1; HPH-3; PHD-1; PHD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001004083
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212035).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004083.1 , NP_001004083.1
RefSeq Size:	2013 bp
RefSeq ORF:	1248 bp
Locus ID:	308457
UniProt ID:	Q6AYU4
Cytogenetics:	1q21
Gene Summary:	HIF-prolyl hydroxylase 1; might be involved in might play a role in hypoxic preconditioning [RGD, Feb 2006]



[View online »](#)