

Product datasheet for **RR212063L3V**

Mthfd1 (NM_022508) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mthfd1 (NM_022508) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Mthfd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022508
ORF Size:	2805 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212063).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022508.1 , NP_071953.1
RefSeq Size:	3332 bp
RefSeq ORF:	2808 bp
Locus ID:	64300
UniProt ID:	P27653
Cytogenetics:	6q24
Gene Summary:	a trifunctional enzyme combining a 10-formyl-H4folate synthetase, a 5,10-methenyl-H4folate cyclohydrolase, and a 5,10-methylene-H4folate dehydrogenase [RGD, Feb 2006]



[View online »](#)