

Product datasheet for **RR212076**

Srsf5 (NM_019257) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Srsf5 (NM_019257) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Srsf5
Synonyms: HRS; Sfrs5; SRp40
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR212076 representing NM_019257
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTGGCTGTCGAGTGTTTCATTGGGAGACTAAATCCCGCAGCGAGGGAGAAAGATGTGGAAAGATTCT
TCAAGGGTTACGGACGGATCAGAGATTTGATTTGAAAAGAGTTTTGGTTTTGTGGAATTTGAGGACCC
AAGGGATCGAGATGACGCTGTTTATGAACCTTGATGGGAAGGAACTTGCAGTGAAAGGGTTACTATTGAA
CATGCTCGGGCTCGGTCTCGAGGTGGAAGAGGTAGAGGACGATACTCCGACCGTTTTAGCAGTCGCAGAC
CTCGAAATGATAGACGAAATGCTCCACCTGTAAGAACAGAAAATCGACTTATAGTTGAGAATTTATCCTC
AAGAGTCAGCTGGCAGGATCTCAAAGATTTTCATGAGACAAGCTGGGGAAGTAACCTTTGCGGATGCACAT
CGACCTAAACTAAATGAAGGGTAGTTGAGTTGCCTCTTATGGTGACTTAAAGAATGCTATTGAGAAAC
TTTCTGGAAGGAAATTAATGGGAGAAAATCAAATTAATTGAAGGCAGCAAAAGACACAGGTCAAGAAG
CAGGTCACGATCTCGACCAGGAGTTCCTCTAGGTCCCGTAGCCGATCCCGTTACGCAGGAGCAAGTCT
TACAGCCGATCAAGGAGCAGGAGCCGGAGCCGGAGCAAGTCCCGTCTGGTAGTAGGTCTCCCGTGCCTG
AGAAGAGCCAGAAACGTGGTCTTCGAGTAGATCTAAATCTCCAGCATCTGTGGATCGCCAGAGGTCCCG
ATCCCGTCCAGGTCCAGATCAGTTGACAGTGGCAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212076 representing NM_019257
Red=Cloning site Green=Tags(s)

MSGCRVFIGRLNPAAREKDVERFFKGYGRIRDIDLKRGFGVFEFEDPRDADDAVYELDGKELCSERVITIE
 HARARSRGGRRGRYSDRFSSRRPRNDRRNAPPVRTENRLIVENLSSRVSWQDLKDFMRQAGEVTFADAH
 RPKLNEGVEFASYGDLKNAIEKLSGKEINGRKIKLIEGSKRHRSSRSRSRTRSSRSRSRSRSRSRSK
 YSRSSRSRSRSKSRGSRSPVPEKSQKRGSSRSKSPASVDRQSRSSRSRSRSVDSGN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_019257
ORF Size: 807 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019257.3 , NP_062130.2
RefSeq Size:	1501 bp
RefSeq ORF:	810 bp
Locus ID:	29667
UniProt ID:	Q09167
Cytogenetics:	6q24
MW:	30.9 kDa
Gene Summary:	may act as a regulator of alternative pre-mRNA splicing; may play a role in cell cycle regulation [RGD, Feb 2006]