

Product datasheet for **RR212088L4V**

Sdpr (Cavin2) (NM_001007712) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sdpr (Cavin2) (NM_001007712) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cavin2
Synonyms:	MGC94054
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001007712
ORF Size:	1251 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212088).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007712.1 , NP_001007713.1
RefSeq Size:	2148 bp
RefSeq ORF:	1254 bp
Locus ID:	316384
UniProt ID:	Q66H98
Cytogenetics:	9q22



[View online »](#)

Gene Summary:

Plays an important role in caveolar biogenesis and morphology. Regulates caveolae morphology by inducing membrane curvature within caveolae. Plays a role in caveola formation in a tissue-specific manner. Required for the formation of caveolae in the lung and fat endothelia but not in the heart endothelia. Negatively regulates the size or stability of CAVIN complexes in the lung endothelial cells (By similarity). May play a role in targeting PRKCA to caveolae (PubMed:9566962).[UniProtKB/Swiss-Prot Function]