

Product datasheet for **RR212114**

Fxr2 (NM_001100647) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fxr2 (NM_001100647) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fxr2
Synonyms:	Fxr2h
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR212114 representing NM_001100647
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCGCCTGGCCTCTGGGGGGATGTGGAGCCGGACTGCCCGTCGAGGTGCGGGGCTCCAACGGGG
 CCTTCTACAAGGCTTTGTTAAGGATGTCCATGAAGACTCTGTCATATCTTCTTTGAAAACAACCTGGCA
 GAGTGAGAGACAGATTCCCTTTGGGGATGTCCGACTACCACTCCAGCTGACTACAATAAGGAGATCACA
 GAAGGAGATGAAGTAGAGGTCTATTCGAGAGCCAATGAACAGGAACCTTGTGGTTGGTGGCTGGCCCGGG
 TGCGGATGATGAAGGGAGATTTCTATGTCATTGAGTATGCTGCCTGTGATGCTACCTACAATGAGATTGT
 CACTCTGGAGCGACTTCGGCCAGTCAATCCCAATCCCCTTGAACCAAAGGCAGCTTCTTCAAGGTTACC
 ATGGCTGTGCCTGAGGATCTTAGAGAGGCTGCTCCAATGAAAATGTCCATAAAGAGTTCAGAAAAGCAC
 TGGGAGCCAACCTGCATCTTCTCAACATCACAGACAGTGAAGCTTTCATTCTGTCAACCACAGAAGCCCC
 TGTGAAGAGGGCATCTTTGCTGGGTGACATGCATTTCCGAAGCCTTCGTACCAAATTCGTGCTCATGTCC
 CGCAATGAAGAAGCTACTAAGCACCTAGAGACAAGTAAGCAATTGGCAGCAGCATTCCAAGAGGAGTTTA
 CAGTTCGAGAGGACCTGATGGGTCTGGCAATTGGGACTCACGGTGCCAACATCCAGCAGGCTCGAAAAGT
 GCCAGGAGTACTGCCATTGAGTTGGGTGAAGAGACCTGCACTTCCGAATCTATGGGGAGACGCCAGAG
 GCATGCAGACAGGCTCGAAGCTACCTTGAATTTCTGAGGACTCTGTACAAGTGGCCAGGAACCTGGTTG
 GCAAAGTGATTGGAAAGAACGGGAAAGTGATTCAGGAGATCGTGGATAAGTCTGGTGTGGTGGCGGTTCCG
 TGTGGAAGGTGATAACGACAAGAAGAACCCAGGGAGGAGGAATGGTTCCCTTCATATTTGTTGGCACC
 CGAGAGAATATCAGCAATGCCAGGCACTGCTAGAGTATCACCTCTCCTACCTGCAGGAGGTTGAACAGC
 TTCGCTTGGAGAGACTGCAAATTGATGAACAGCTTCGGCAGATTGGGCTAGGCTTTCGCCCTCTGGGAG
 TGGACGGGGTGGCAGTGGTGGTGGCAGTGACAAGGCGGGCTATACCACAGACGAGAGCTCCTCCTCTTCC
 CTGCACACGACCAGAACCTATGGGGCAGCTATGGGGCCGGGGCCGGGGCCGGAGGACTGGTGGTCTG
 CCTATGGCCCCAGCTCAGACCCCTCCACAGCTTCTGAGACAGAGTCAGAAAAGAGAGAGGAATCTAATCG
 AGCAGGACCTGGTGACAGAGATCCCCAAGCCGTGGAGAAGAAAGCAGGAGGCGGCAATTGGGGCCGG
 GGAAGGGGACCCCACTGTTCCCCGGCTACCTCAAGATACAACCTCATCATCTATTAGCTCAGTGCTAA
 AGGATCCAGACAGCAACCCTTACAGCCTATTGGACACATCTGAACCAGAGCCTCCAGTTGATTGAGAACC
 TGGTGAACCTCTCCAGCAAGTGCCCGTCGTCGCCGTTCCCGCCCGCTGCTACTGATGAAGACAGGACT
 GTCATGGACGGAGCCTGGAGTCTGACGGTCTAACGTGACAGAGAATGGCCTGGAAGATGAATCAAGAC
 CCCAGAGACGTAATCGCAGCCCGCCGCGGTAATCGTGGTAACCGGACTGATGGGTCCATCAGCGGAGA
 CCGCAACCAGTACTGTGGCTGACTACATCTCACGGCAGAGTCTCAGAGCCGCCAGCGGCCACCTCTG
 GAGCGCACAAAACCTCAGAGGATTCTTTTCAGGACAGAAGGGTACTCTGTCAGCAGGCTTCCCAAGG
 GCCCTCAGAGAATGGGGAGCTCTGCCCCATTGGAGTTGGGTAGTTTGGTGAATGGGGTTTCA

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR212114 representing NM_001100647
Red=Cloning site Green=Tags(s)

```

MGGLASGGDVEPGLPVEVRGNSNGAFYKGFVKDVHEDSVTIF FENNWQSERQIPFGDVRLPPPADYNKEIT
EGDEVEVYSRANEQEPCGWLLARVRMMKGDFFYVIEAACDATYNEIVTLERLRPVNPNPLATKGSFFKVT
MAVPEDLREACSNENVHKEFKKALGANCIFLNITDSEF ILSTTEAPVKRASLLGDMHFRSLRTRKLLMS
RNEEATKHLETSKQLAAAFQEEFTVREDLMGLAIGTHGANIQARKVPGVTAIELGEETCTFRIYGETPE
ACRQARSYLEFSEDSVQVPRNLVGKVIKNGKVIQEIVDKSGVVRVVRVEGDNDKKNPREEGMVPIFVGT
RENISNAQALLEYHLSYLQEVEQLRLERLQIDEQLRQIGLFRPPGSGRGGSGGSDKAGYTTDESSSS
LHTTRTYGGSYGRGRGRRTGGPAYGPSSTPASETSEKREESNRAGPGDRDPPSRGEESSRRRPIGGR
GRGPPVPRPTRSYNSSISSVLKDPDSNPYSLLDTSEPEPPVDESEPPPPASARRRRRRRRRTDEDRT
VMDGGLESDPNVTENGLEDESRPQRNRNRNRNRNRNRDGSISGDRQPVTVADYISRAESQSRQRPPL
ERTKPSSEDSLGGQKGDVSRLPKGPSENGELSAPLELGSLVNGVS
  
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001100647

ORF Size: 2025 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001100647.1](#), [NP_001094117.1](#)

RefSeq Size: 2671 bp


RefSeq ORF: 2028 bp

Locus ID: 287433

Cytogenetics: 10q24

MW: 74.4 kDa

Product images:



Circular map for RR212114