

Product datasheet for **RR212118L3V**

Pcsk9 (NM_199253) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pcsk9 (NM_199253) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pcsk9
Synonyms:	NARC-1; Narc1; PC9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199253
ORF Size:	2073 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212118).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199253.2
RefSeq Size:	3350 bp
RefSeq ORF:	2076 bp
Locus ID:	298296
UniProt ID:	P59996
Cytogenetics:	5q34
Gene Summary:	a soluble zymogen; involved in hepatic growth and differentiation [RGD, Feb 2006]



[View online »](#)