

Product datasheet for **RR212187**

Ccnd1 (NM_171992) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccnd1 (NM_171992) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ccnd1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR212187 representing NM_171992
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAACACCAGCTCCTGTGCTGCGAAGTGGAGACCATCCGCCGCGCGTACCCTGACACCAATCTCCTCA
ACGACCGGTGCTGCGAGCCATGCTTAAGACTGAGGAGACCTGCGGCCCTCCGTTTCTTACTTCAAGTG
CGTGCAGAGGGAGATTGTGCCATCCATGCGGAAAATCGTGCCACCTGGATGCTAGAGGTCTGCGAGGAG
CAGAAGTGCGAAGAGGAGGTCTTCCCGCTGGCCATGAACACTACCTGGACCGTTTCTTGTCTCTGGAGCCCC
TGAAGAAGAGCCGCTGCAGCTTCTGGGGCCACCTGCATGTTCTGGCCCTAAGATGAAGGAGACCAT
TCCCCTGACTGCCGAGAAGTTGTGCATCTACACTGACAACCTATCCGCCCGAGGAGTTGCTGCAAATG
GAACTGCTTCTGGTGAACAAGCTCAAATGGAACCTGGCCGCCATGACTCCCCACGATTTTCATCGAACACT
TCCTCTCCAAAATGCCAGAGGCGGATGAGAACAAGCAGATCATCCGAAACATGCACAGACCTTTGTGGC
CCTCTGTGCCACAGATGTGAAGTTCATTTCCAACCCACCCTCCATGGTGGCTGCTGGGAGTGTGGTGGCC
GCGATGCAGGGCCTGAACCTGGGAGCCCCAACAACTTCTCTCCTGCTACCGCACAAACGCACTTTCTTT
CCAGATCATCAAGTGTGACCCGGACTGCCTCCGTGCCTGCCAGGAACAGATTGAAGCCCTTCTGGAGTC
AAGCCTGCGCCAGGCCAGCAGAATCGATCCCAAGGCCACCGAGGAGGAAGGGGAAGTGGAGGAAGAA
GCTGGTCTGGCTTGCACACCCACCGACGTCAGATGTGGACATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR212187 representing NM_171992
Red=Cloning site Green=Tags(s)

```
MEHQLLCCEVETIRRAYPDTNLLNDRVLRAMLKTEETCAPSVSYFKCVQREIVPSMRKIVATWMLEVCEE
QKCEEEVFPLAMNYLDRFLSLEPLKKSRLQLLGATCMFVASKMKETIPLTAEKLCIYTDNSIRPEELLQM
ELLLVNKLKWNLAAMTPHDFIEHFLSKMPEADENKQIIRKHAQTFVALCATDVKFI SNPPSMVAAGSVVA
AMQGLNLGSPNNFLSCYRTHFLSRVIKCDPDCLRACQEQIEALLESSLRQAQQNIDPKATEEEGEVEEEE
AGLACTPTDVRDVDI
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1865_g02.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_171992


ORF Size: 885 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_171992.4</u> , <u>NP_741989.3</u>
RefSeq Size:	3809 bp
RefSeq ORF:	888 bp
Locus ID:	58919
UniProt ID:	<u>P39948</u>
Cytogenetics:	1q42
MW:	33.9 kDa
Gene Summary:	regulatory subunit of CDK4 or CDK6; involved in regulation of the G1/S transition of the cell cycle [RGD, Feb 2006]

Product images:



Circular map for RR212187