

Product datasheet for **RR212203L3V**

Epha7 (NM_134331) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Epha7 (NM_134331) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Epha7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134331
ORF Size:	2994 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212203).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134331.1 , NP_599158.1
RefSeq Size:	3208 bp
RefSeq ORF:	2997 bp
Locus ID:	171287
UniProt ID:	P54759
Cytogenetics:	5q21
Gene Summary:	receptor for members of the ephrin-A subfamily of receptor tyrosine kinases; binds with ephrin-A1, -A2, -A3, -A4 and -A5, may mediate developmental effects in the nervous system [RGD, Feb 2006]



[View online »](#)