

Product datasheet for **RR212216L4V**

Ptger2 (NM_031088) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ptger2 (NM_031088) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ptger2
Synonyms:	EP2; Ptger-ep2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031088
ORF Size:	1071 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212216).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031088.1 , NP_112350.1
RefSeq Size:	1360 bp
RefSeq ORF:	1074 bp
Locus ID:	81752
UniProt ID:	Q62928
Cytogenetics:	15p14
Gene Summary:	binds prostaglandin E2 and induces increased cAMP; may play a role in skeletal growth [RGD, Feb 2006]



[View online »](#)