

Product datasheet for **RR212226L3V**

Rfc1 (NM_053547) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rfc1 (NM_053547) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rfc1
Synonyms:	Recc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053547
ORF Size:	3393 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR212226).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053547.1 , NP_445999.1
RefSeq Size:	4741 bp
RefSeq ORF:	3396 bp
Locus ID:	89809
Cytogenetics:	14p11
Gene Summary:	transcriptional repressor that regulates transcription of the vasoactive intestinal peptide receptor gene [RGD, Feb 2006]



[View online »](#)